Microsart® e.jet

Laboratory Transfer Pump with Quick Connection for Microbiological Analysis

Microbiological Analysis

Typical Set-up: Microsart® e.jet connected to a 3-branch manifold with Microsart® Funnels 100

New pump for vacuum filtration in Microbiological Analysis

The Microsart® e.jet is a laboratory vacuum pump which creates vacuum and concomitantly transfers the filtered liquid to waste. The Microsart® e.jet is ideal for sample preparation in microbiological analysis achieving a trans membrane pressure of 600 mbar and a flow rate of 4.0 NI/min (4.0 Normliters water displacement by air in one minute). Constant flow rates and a defined maximum vacuum guarantees a smooth and reliable filtration.

Reduced operating complexity

Until now vacuum equipment for the Membrane Filtration Method consists of numerous parts including connectors, tubes, vacuum containers, protection filter, Woulff's bottle and a vacuum pump. After several samples the vacuum must be broken to empty the filtrate collection container. The complete traditional equipment requires far more laboratory space and is time consuming to operate and maintain. The new Microsart® e.jet greatly reduces operating complexity.

Smart design fits any laboratory environment

The Microsart® e.jet pump is an ideal accessory for 3-branch and 1-branch manifolds. Compared to traditional equipment, the

Microsart® e.jet and stainless steel manifold require only 30% of the average space, which translates to less congestion when working in laminar flow cabinets.

Safe operation | Maintenance-free

Traditional vacuum pumps often lose their efficiency and capability to generate sufficient vacuum, when liquid is drawn into the pump head. The Microsart® e.jet is designed to pump both gas and liquids, meaning no loss of efficiency or malfunction from water drawn into the pump head.

Building-up the vacuum filtration system is easy and fast thanks to the innovative Quick Connections. The Microsart® e.jet Transfer Pump is equipped with Quick Connection Nipples assembled to Quick Connection Couplings on hose nipples for DN 10 tubings. Simply push-to-connect for assembling and pull-to-disassembling the whole system within seconds. The Quick Connections are non-shut-off.

Technical Specification

Flow rate	> 4.0 NI/min
Max. Vacuum	0.4 bar
Max. Pressure	1.0 bar
Materials (contact with filtrate)	PTFE, ETFE, Polypropylene, EPDM, POM, PSU
Mains	100-240 V 47-63 Hz
Weight Pump	1,425.3 g
Weight Power supply	242.6 g
Dimensions (W x L x H)	12 x 17 x 19 cm
Max. ambient temp.	+5...+40°C
Max. temp. of liquid	+5...+80°C
Max. viscosity	< 150 cSt*
Protection type	IP 64
Protection class	III
Inlet Outlet	Quick Connection on hose nipples for DN 10 tubings

Ordering Information

Order Number	Description
166MP-4	Microsart® e.jet Transfer Pump

Spare Parts

Order Number	Description
1EP---0003	Pump head complete for 166MP-4
1EE---0012	Power supply complete for 166MP-4
1EAS--0027	2 Quick Connection Couplings (PSU) and 2 Nipples (POM)
1ZE---0053	Foot switch

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178
UK +44.1372.737159
France +33.442.845600
Italy +39.055.63.40.41
Spain +34.90.2110935
Japan +81.3.3740.5407

Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH. | W
Publication No.: SM-2008-a11094
Order No.: 85034-539-02
Ver. 09 | 2011

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