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Directions for Use

16824

3-Branch Filtration System with 100 ml Funnels

16828

3-Branch Filtration System with 500 ml Funnels

16831

6-Branch Filtration System with 500 ml Funnels

16832

6-Branch Filtration System with 100 ml Funnels

Contents

1. Description of the Filtration Systems
2. Specifications
3. Equipment Supplied
4. Cleaning and Care
of a Multibranch Filtration System
5. Microbiological Tests
6. Analytical Procedures
7. Accessories
8. Spare Parts

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1. Description of the Filtration Systems

Each multibranch system consists of a manifold with 3 or 6 filter holder support bases welded in place, each complete with a silicone gasket below the stainless steel filter support frit, and a stainless steel tap for independent vacuum control for each funnel. All funnels have a quick-lock clamp, for simple and secure mounting on their respective support bases, and a lid with an inner silicone sealing ring.

The stainless steel funnels also have convenient level marks to simplify filtration of premeasured quantities. They are located on the inside opposite to the clamp at levels of approx. 100, 200, 300, 400 and 500 ml for 500 ml funnels and at levels of approx. 50 and 100 ml for 100 ml funnels. Please check the real volume of your sample with a suitable volume measure device.

Advantages:

- Timesaving for routine analyses: 3 or 6 samples, resp., can be filtered simultaneously
- Durable
- Practical and easy to use
- Stainless steel taps permit independent control of vacuum supplied to each funnel.

Applications:

- Microbiological quality control, e.g., detection of microbes in water and other beverages
- Analyses

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Vacuum Connection

The manifolds have a hose barb with an 10 mm outer diameter to attach a vacuum hose, such as the Sartorius Stedim Biotech rubber vacuum hose 16623 for connection to a vacuum source.

The vacuum source can be an installed vacuum supply pipeline or

- a) a water jet pump (e.g., 16611, with an R 3/4" female thread)
- b) an electric pump (e.g., laboratory vacuum pumps, 16612 or 16692, 220 V 50 Hz, also available for 110 V 60 Hz under order numbers 16615 and 16695, resp.) In this case, a vacuum resistant flask must be interconnected between the pump and manifold to collect the filtrate.

Prefiltration

Samples with a heavy "dirt load" or a high content of insoluble solids must be prefiltered for bacteriological analysis. Stainless steel prefilter attachment 16807 lets you do prefiltration and actual membrane filtration in one easy step.

The attachment with a bacteriological prefilter in place, 12301-050 (8.0 µm pore size), is simply mounted between the funnel and the base. Water-insoluble, coarse particulates in the sample are retained on the prefilter, allowing the bacteria-retentive membrane filter on the underlying support to trap bacteria without an interfering buildup of particulate matter that may hide or disturb the growth of colonies.

The filter support can be flamed as is possible with the filter support bases on the vacuum manifolds. All filter holders are also autoclavable at 121°C.

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Stainless steel

Lid sealing

16824 and 16832:
silicone sealing ring, 77.2/85.8 mm
16828 and 16831:
silicone sealing ring, 122/131 mm

Funnels

Complete with stainless steel clamp.
16824: 3 funnels,
16832: 6 funnels, each with a 100 ml
capacity
16828: 3 funnels,
16831: 6 funnels, each with a 500 ml
capacity

Filter supports

Stainless steel frit, 50 mm

Gaskets

Silicone gasket

Manifold base

Stainless steel. Equipped with 3 (16824
and 16828) or 6 (16831 and 16832) filter
support bases with taps, and a hose barb for
10 mm nominal diameter hoses

Suitable filter diameter

Membrane filters: 50 mm (or 47 mm after
exchanging the supplied stainless steel frit
for the 47 mm accessory frit 6980103)

Operating pressure

Vacuum only

Sterilization

Autoclaving up to 134°C or dry heat at
180°C without filters in place; sanitization
by flaming

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3. Equipment Supplied

Each multibranch manifold is supplied completely assembled. The taps and clamps are lubricated for your convenience.

4. Cleaning and Care of a Multibranch Filtration System

4.1 The multibranch filtration system must be cleaned and dried before and after use. Flush it with hot water and rinse with distilled water. To remove stubborn stains, dirt, and other impurities, you may additionally use a commercially available laboratory detergent (for metal, glass, and plastic) and soft brushes.

4.2 Care of the taps:

4.2.1 The taps come pregreased (we use high-grade, heavy grease). Therefore, upkeep involving disassembly, cleaning, drying, and lubrication should not be necessary unless they later become difficult to open and close. After cleaning and greasing the taps, make sure to install each one in its original tapered seating.

4.2.2 Disassembly of the taps:

- a) Use flat nose pliers or 3.2 mm Bensing gripping pliers to remove the lock washer.
- b) Remove the leaf spring and the washer by hand.
- c) Pull out the tap to remove.

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4.2.3 Installation of a tap:

Clean, dry and regrease the tap. Make absolutely sure that you install it in its original stem guide (tapered seating) since the valve stem guide has been precisely machined to the tap's individual specifications to make a leaktight fit. If you randomly install the taps, they may leak or not remain in the desired position.

4.3 Care of the clamps:

4.3.1 The clamps tend to widen after frequent use and must be retightened to prevent leakage.

4.3.2 To retighten a clamp, grasp it with flat nose pliers, applying slight pressure to narrow its opening.

4.3.3 The moving parts of the clamp are exposed to friction caused by opening and closing. Always keep them lightly greased.

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5. Microbiological Tests

5.1 Clean your laboratory benchtop and disinfect it with alcohol or a commercially available disinfectant.

5.2 Line up your sterile culture media in a row, such as nutrient pads in petri dishes. Place the packaged sterile membrane filters within easy reach.

5.3 Fill a glass beaker with a slight amount of alcohol and place forceps inside. Position it along with a Bunsen burner next to the multibranch filtration system (Fig.1).

5.4 Switch on your vacuum source (pump or water jet pump) and remove a stainless steel funnel from the manifold (Fig. 2).

5.5 Flame the filter support base and the metal frit by opening the tap to draw the flame through them. Close the tap (Fig.3).

5.6 Now pick up the stainless steel funnel and flame it on the bottom (Fig.4).

5.7 Set the stainless steel funnel on the filter support base and flame the inside from the bottom to the top, moving the flame in an upward spiral (Fig.5).

5.8 Flame the lid (Fig.6) and position it on top of the funnel. Now the filter holder is sterile. Repeat this procedure to prepare the other filter holders.

5.9 Using Sartorius Stedim Biotech blunt-tipped forceps, which you have sterilized by brief flaming, pick up a sterile membrane filter; use your free hand to lift the funnel together with the lid, and center the membrane filter on the frit (Fig.7).

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5.10 Now reposition the funnel and secure in place with the clamp (Fig. 8). Follow this procedure to prepare the remaining filter holders of your 3 or 6-branch filtration system.

5.11 After pouring the samples in each of the stainless steel funnels, open the taps to filter.

5.12 At the end of the filtration run, close the taps. Use flamed forceps to remove the filters from the filter holders, and place each one with the grid faceup on a moistened nutrient pad (Fig. 9).

5.13 Next, incubate the petri dishes in an incubator as directed.

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6. Analytical Procedures

6.1 If necessary, flush the entire filtration assembly free of particles. In addition, rinse your forceps and membrane filters to wash off any particles.

6.2 Unmount a funnel, and use Sartorius Stedim Biotech blunt-tipped forceps to center a membrane filter on the frit.

6.3 Now reposition the funnel and secure it in place with the clamp. Prepare the remaining filter holders of the 3- or 6-branch manifold in the same manner.

6.4 Pour the samples into the respective funnels.

6.5 Switch on the vacuum source and open the taps to filter the samples.

6.6 At the end of the filtration run, close the taps, remove the stainless steel funnels and then the membrane filters to perform your particular analytical procedure (microscopic, gravimetric, chemical or other analysis).

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7. Accessories


16807	Stainless steel prefilter attachment
16623	Rubber vacuum hose (1 m)
16611	Water jet pump with R ¾" female thread
16612	Laboratory vacuum pump, max. 26 Lpm, 220 V, 50 Hz
16615	Laboratory vacuum pump, max. 26 Lpm, 110 V, 60 Hz
16692	Vacuum pump, max. 20 Lpm, 220 V, 50 Hz
16695	Vacuum pump, max. 20 Lpm, 110 V, 60 Hz
18119	Incubator, 100-240 V, 50/60Hz, 18 Litre
16625	Stainless steel forceps for membrane filters
16610	Woulff's bottle, 500 ml
16672	Vacuum flask, 2 l
16672-----1	Vacuum flask, 5 l
17173	Stopper for 16672
17649	Colony Counter
6980104	PTFE gasket (2/pack)

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8. Spare Parts

6981002	Stainless steel funnel, 500 ml, with clamp
6981065	Stainless steel funnel, 100 ml, with clamp
6981004	Stainless steel funnel, 40 ml, with clamp
6980103	Stainless steel frit (47 mm)
6981003 6981064	Silicone sealing ring, 122/131 mm Silicone sealing ring, 77.2/85.8 mm
6980102	Stainless steel frit (50 mm)
6980124	Silicone gasket (2/pack)
6980267	Stainless steel tap
6980271	Silicone O-ring, 16x2 mm
6980268	Knurled nut
6981001	lid for funnel, 500 ml
6981063	lid for funnel, 100 ml
6982119	Spare part Kit for installing the tap (consists of 1 leaf spring, 1 washer, 1 lock washer)
6980270	Stainless steel cap
6980264	Klammer

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
3-Branch Filtration System with 500 ml Funnels (left), with 100 ml Funnels (center) and with 40 ml Funnels (right)

1

2

3

4


5

6

7

8

9

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.

All rights reserved. No part
of this publication may
be reprinted or translated in
any form or by any means
without the prior written
permission of Sartorius Stedim
Biotech GmbH.

The status of the information,
specifications and illustrations
in this manual is indicated
by the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications and
design of the equipment
without notice.

Status:
September 2014,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany