



Preventive Maintenance BIOSTAT® Bioreactor|Fermenter Order No.: S87**SM

Service



The objective of preventive maintenance is a continuous, precise operation of your equipment, to secure your investments for the long term and to guarantee correct and proper function.

Preventive maintenance will be performed in accordance with Sartorius' fixed standards and in between the contractually stipulated time intervals between the customer and Sartorius.


Specifications:

- Inspection visit at fixed intervals or as needed
- Visual inspection
- Maintenance check according to maintenance checklist
 - Pipe system
 - Valve
 - Vessel
 - Pump
 - Motor
- According to maintenance checklist, exchange of
 - Gaskets and o-rings
 - Valve membranes
 - Flexible hose connections
 - Blower filter
- Start-up test
- Leakage test with pressure (only for in situ sterilizable fermenters)
- Sterilization (only for in situ sterilizable fermenters)

Prerequisites required by the customer:

- Unobstructed access to the system
- When ordering, a responsible person has to be appointed
- Access by automobile is possible | parking lot is available

Repairs, spare parts as well as shipping costs, travel costs and expenses are not covered in scope of services and will be invoiced separately.



Device-specific Order Numbers	Order No.
Preventive maintenance BIOSTAT® A plus	S8710SM
Preventive maintenance BIOSTAT® B, B plus	S8713SM
Preventive maintenance BIOSTAT® B plus Twin	S8759SM
Preventive maintenance BIOSTAT® B-DCU II Single	S8761SM
Preventive maintenance BIOSTAT® B-DCU II Twin	S8762SM
Preventive maintenance BIOSTAT® B-DCU II Triple	S8763SM
Preventive maintenance BIOSTAT® B-DCU II Quad	S8764SM
Preventive maintenance BIOSTAT® B-DCU II 5-fold	S8794SM
Preventive maintenance BIOSTAT® B-DCU II 6-fold	S8795SM
Preventive maintenance BIOSTAT® C plus	S8760SM
Preventive maintenance BIOSTAT® D	S8715SM
Preventive maintenance BIOSTAT® D-DCU	S8784SM
Preventive maintenance BIOSTAT® PBR2S	S8785SM
Preventive maintenance BIOSTAT® Q plus 3	S8769SM
Preventive maintenance BIOSTAT® Q plus 6	S8770SM
Preventive maintenance BIOSTAT® Q plus 9	S8771SM
Preventive maintenance BIOSTAT® Q plus 12	S8772SM
Preventive maintenance BIOSTAT® Cultibag STR 500L	S8787SM
Preventive maintenance BIOSTAT® Cultibag STR 1000L	S8788SM
Preventive maintenance BIOSTAT® Cultibag STR plus 50 200	S8765SM
Preventive maintenance BIOSTAT® Cultibag STR plus 50 200 Twin	S8766SM

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

Specifications subject to change without notice.

www.sartorius-stedim.com

Date: August 2012